

Epo-Grip



Fin Magic is a time saving reinforced translucent fin reproduction adhesive designed to produce strong, flexible fins for the taxidermist without the mess. Fin Magic is a ready to use product that can be easily tinted, providing application flexibility not previously found in other methods fin reproduction.

Directions:

Coat fin replica mold with an appropriate release material.

Mix equal portions of Fin Magic sides A & B until a uniform consistency is obtained. Apply Fin Magic to each side of fin mold and clamp sides securely until cure occurs in approximately 12 to 15 minutes. Cure time will vary with the thickness of the fin being made. Thicker fins will cure faster than thin reproductions. After cure separate mold halves, remove the Fin Magic reproduction. Trim to desired size and shape.

Coloring:

Fin Magic can be easily be pre-colored by adding shaved artist pastel chalk to the mixture and blending to desired color.